INSTRUCTIONS FOR USE

NeoSEALER® EZFlo

ENDODONTIC SEALER

COMPOSITION & DESCRIPTION

NeoSEALER EZFlo is a premixed bioactive bioceramic root canal sealer formulated with tricalcium and dicalcium silicate cements for permanent root canal obturation. NeoSEALER EZFlo sets in vivo in the presence of moisture provided by the surrounding tissues. Optimized for easier dispensing through a 29ga tip, NeoSEALER EZFlo is designed to support minimally invasive endodontic procedures which respects the natural anatomy of the canal.

CLINICAL BENEFITS/PERFORMANCE CHARACTERISTICS

- Bioactive bioceramic
- Does not discolor teeth
- Radiopaque
- Resin-free
- Easier dispensing using a 29ga Tip
- Optimized for both single cone and warm obturation techniques

INTENDED USE, INDICATIONS, TARGET GROUP AND INTENDED POPULATION

Avalon Biomed NeoSEALER EZFlo is designed for obturation and sealing of root canals following removal of vital pulp tissue, removal of infected/necrotic pulp and placement of intracanal dressings, or in previously treated teeth requiring endodontic retreatment. Avalon Biomed NeoSEALER EZFlo is for pediatric and adult patients (2 years of age or older) undergoing root canal therapy on permanent teeth with closed apices, or on primary teeth with no permanent successor. Avalon Biomed NeoSEALER EZFlo is for professional (Dentists/Endodontists) use only, to be used in clinical environments (hospitals, clinics & dental offices).

SUBSTANCE EXPOSURE:

The premixed root and pulp material is composed of modified ceramic powder mixed with an organic liquid to a low viscosity paste and may come into incidental contact with mucosal membranes. The main constituents of Avalon Biomed NeoSEALER EZFlo are as follows. Refer to the product Safety Data Sheet for additional information.

Raw Materials	CAS. Number	EC Number	Content (%)
Ta ₂ O ₅ -Tantalum pentoxide	1314-61-0	215-238-2	<50
3CaO·SiO ₂ -Tricalcium Silicate	12168-85-3	235-336-9	<25
CaAl ₂ O ₄ -Calcium Aluminate	65997-16-2	266-045-5	<25
2CaO·SiO ₂ -Dicalcium Silicate	10034-77-2	233-107-8	<10
CaAl ₂ O ₇ -Calcium Dialuminate	12004-88-5	234-463-7	<15
Al ₂ Si ₂ O ₇ (OH) ₂ ·2H ₂ O-Aluminosilicate	1332-58-7	310-194-1	<10

CONTRAINDICATIONS

- Hypersensitivity against caustic (high pH) solutions.
- Do not use for primary tooth pulpectomy (root canal filling) unless the permanent successor tooth is absent.

WARNINGS

NeoSEALER EZFlo is caustic, as are all tricalcium silicates.

PRECAUTIONS

- NeoSEALER EZFlo is a non-cytotoxic, biocompatible paste. Resorption in the periapical tissue may occur over years. Caution should be exercised to NOT overfill the root canal forcing a large amount of sealer beyond the apex of the root. If a large amount of sealer is overfilled into the mandibular canal (inferior alveolar canal), the clinician should assess if the overfilled material should be removed according to state-of-the-art policy.
- Minimize over extension of the material beyond the apex.
- Avoid contact of unset sealer with skin or oral mucosa. After incidental contact, wash and rinse with water.
- Wear suitable gloves and protective glasses during use.
- NeoSEALER EZFlo MUST BE KEPT WELL SEALED. Immediately recap after each use.
- To protect against moisture intrusion, store NeoSEALER EZFlo in its protective foil pouch.
- Always use a new tip when applying NeoSEALER EZFlo.
- Avoid contacting the syringe to a contaminated surface.
- Cover the syringe body with a disposable protective sleeve if used intraorally, to minimize contamination of the syringe.
- NeoSEALER EZFlo and tips are provided in non-sterile packaging. These products cannot be sterilized. Clinicians should follow their established protocols for cleaning and disinfection of the NeoSEALER EZFlo syringe between uses.

See: <https://www.cdc.gov/infection-control/media/pdfs/Guideline-Disinfection-H.pdf>

- Setting of tricalcium silicates may be inhibited in acidic environments such as infected sites.

ADVERSE REACTIONS

Reversible acute inflammation of the oral mucosa may occur if contacted with the unset paste.

MALFUNCTIONS & SERIOUS INCIDENTS:

If any malfunction or changes in performance that may affect safety occur, the clinician should assess and take appropriate action deemed necessary for the patient if affected by the malfunction or changes in performance. Clinician should report malfunctions and/or changes in performance that may affect safety to NuSmile, Ltd. Any serious incident that has occurred in relation to the device should be reported to Avalon Biomed via email (info@avalonbiomed.com) or phone (+1-941-896-9948) and the Competent Authority of the Member State in which the user and/or patient is established.

SYMBOLS USED ON LABELING:

	Manufacturer		Keep Dry		Medical Device		Authorized Representative in the Switzerland Community
	Unique Device Identifier		Lot Number		Non-Sterile		Authorized Representative in the United Kingdom Community
	Prescription Only		Catalog Number		Consult Instructions For Use		Authorized Representative in the European Community
	Irritant		Expiration Date		Caution	 Mfg. by NuSmile, Ltd 3315 West 12th St Houston, TX 77008 USA +1.713.861.0033	

Non-Staining
BIOACTIVE
Bioceramic

NeoSEALER® EZFlo

INTERACTIONS WITH OTHER DENTAL MATERIALS

None known.

STORAGE

Store at room temperature. (DO NOT REFRIGERATE). To prevent hardening of the NeoSEALER EZFlo, immediately recap after each use. Store the syringe in the protective foil pouch provided.

CLINICAL DIRECTIONS FOR USE:

- Excavate carious tooth structure using a round bur in a handpiece at low speed or use hand instruments.
- Debride, clean, shape, and disinfect the root canal system using intra-canal instruments under rubber dam isolation.
- After preparation and irrigation of the root canal, including removal of smear layer, dry the canal using paper points to remove pooled irrigants or body fluids. Moisture is required for setting.
- NeoSEALER EZFlo may be used with or without the tip, but the tip is recommended.

USE WITHOUT the Tip:

- For obturation techniques where most of the canal is obturated by gutta percha points, first apply a light coating of NeoSEALER EZFlo to the canal walls using a device of your choice.
- Second, coat the disinfected and dried gutta percha points with NeoSEALER EZFlo material before insertion.

USING the Tip:

- Attach the tip to the NeoSEALER EZFlo syringe.
- Insert the tip into the canal to the desired depth, taking care to not insert the end past the middle third of the root canal.
- Do not force the tip into the canal such that it seals against the canal walls. Doing so may force the sealer to be expressed beyond the apex.
- Apply gentle, steady pressure to the syringe plunger with a constant outward movement from the canal as the material is expressed.

NOTE: NeoSEALER EZFlo is optimized for use with a 29ga Tip. Other compatible dispensing tips may be used at the clinician's discretion. Smaller diameter tips (higher gauge sizes) may increase dispensing resistance.

- Slowly, insert a single or multiple gutta percha cones into the canal.
- Melt and cut off excess gutta percha at the orifice following the manufacturer's directions.
- Remove excess sealer with a moist cotton pellet.
- Confirm placement of the material in the complete root canal system and check for any over-extrusion with a radiograph.

OTHER APPLICATIONS:

NeoSEALER EZFlo may be used with NeoPUTTY® for obturation of primary teeth with no successor, or otherwise unsalvageable teeth.

NOTE: For removal of root canal filling, use conventional instrumentation to remove set sealer, including hand/rotary files and ultrasonic instrument tips, as with any other sealer. A bioceramic sealer solvent may also be used to assist in removal.

SUMMARY OF SAFETY CLINICAL PERFORMANCE:

<https://ec.europa.eu/tools/eudamed/#/screen/home>. This webpage will be replaced with the direct link to NuSmile's Summary of Safety Clinical Performance (SSCP) when available on EUDAMED's portal by the European Union Regulatory Authority. NuSmile's Summary of Safety Clinical Performance may be provided upon request.

MISCELLANEOUS:

IFU-128 Avalon Biomed NeoSEALER EZFlo Instructions for Use is available electronically at www.avalonbiomed.com. Physical copies shall be provided upon request, free of charge.

ADA 57, ISO 6876

Test	Acceptance Criteria	NeoSEALER EZFLO Results
Working time at room temperature:	None	~35 minutes
Setting time at 37°C, in vivo (or moist environment):	Material must set	Meets ISO 6876 standard
Flow:	>17 mm	Meets ISO 6876 standard
Film Thickness:	<50 µm	<50 µm
Solubility:	<3% by weight	<3% by weight
Dimensional stability:	<1% shrinkage & <0.1% expansion	~0.04% expansion
Radiopacity:	>3 mm equivalent of Al	5.5 mm of Al

IFU-128 Rev 1 Avalon Biomed NeoSEALER EZFlo
English
Effective Date: 07/14/2026

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